Exceptional ergonomics by design

DP-50 has been designed following a careful study of the changing needs of clinical practices. A smart new shape, enhanced mobility and more convenient operation make it a system that is well-suited for clinical settings.

- Sleek, streamlined, compact shape
- User friendly keyboard and controls
- Lightweight system for enhanced mobility
- High resolution, wide-angle 15" LCD with tilt functionality for better viewing



DP-50

Digital Ultrasonic Diagnostic Imaging System


As a stepping stone into the world of ultrasound, B/W ultrasound today need to bring more advances to general practices to improve their patient's benefits. Being a world-class medical equipment provider, Mindray's deep understanding of needs from customers makes our system DP-50 more suitable for your daily diagnosis in less time, while maintain your budget needs.













Diagnostic confidence with advanced imaging technologies

The DP-50 incorporates advanced technologies that deliver deeper penetration and higher resolution. This results in a more confident diagnosis across a range of clinical applications.

iTouch Auto Optimization


One touch image optimization, quickly optimizes the uniformity and dynamic range, very convenient and efficient






iClear Speckle Suppression

Scale the image size to intelligently detect structures of different sizes, improve accuracy



Phase Shift THI

Enhance contrast resolution and improve image quality by separating basal frequency and harmonic content



Efficient your workflow intelligently

DP-50's concise graphical user interface is well organized as control area, clip board area and imaging area. This helps users to remain focused and to leverage the full diagnostic potential of the system all the times.





Simplified control area, whatever operation you do, comment added, measure, image control or setting parameters, just select then get result.

2. Clip Board Area



Intelligent clip board, Cineloop or frame can be easily previewed, played, sent or deleted.

3. Imaging Area



Large imaging area offer more information to make accurate diagnosis

Quick Save

One push image transfer to local or directly to USB

Onboard reporting

Standard templates for various applications

DICOM 3.0 Compliant

Structured reporting Worklist

DICOM print

Storage



DP-50

Digital Ultrasonic Diagnostic Imaging System

The materials and information provided via this brochure are distributed internationally, but not all the products or services found on this brochure are necessarily available in your country or area. Please contact local office or agents to make sure of the availability of the relevant products or services.

Address: Rua Joaquim Floriano, 488, floor Intermediário - cj 02, Itaim Bibi, São Paulo - SP, CEP 04534-002, Brazil Tel: (55-11) 3078-8026

Address: 4250 Kingsway, Suite 206, Burnaby, BC, Canada V5H 4T7 Tel: (1-604) 451-1199

EgyptAddress: 24th floor,34 EL-Obour Buildings, Salah Salem Street, Helioplis, Cairo, Egypt Tel: (20-2) 2260-9728

Address: Europarc Créteil Cedex 123, Chemin des Bassins 94000 Créteil France Tel: (33-1) 4513-9145

Germany Address: Zwischen den Bächen 4, 64625 Bensheim Tel: (49-6251) 17524-0

India
Address: Unit No. 401/402, NDM-1, Netaji Subhash PlaceWazirpur District Centre, New Delhi 110 034 India Tel: (91-11) 4923-0000

Indonesia Address: Menara Dea 1, Suite 301, Kav E4 No.3/1, Mega Kuningan, Jakarta Selatan 12950, Indonesia Tel: (62-21) 576-2650

Italy Address: Via Girardo Patecchio, 4 - 20141 Milano, Italy Tel: (39-02) 57402-4444

Address: Insurgentes Sur # 1787 Piso 2, Col. Guadalupe Inn, Mexico D.F. 01020 Tel: (52-55) 5661-9450

Netherlands Address: Saturnusstraat 17, 2132HB Hoofddorp, The Netherlands Tel: (31-23) 574-9656

Address: Kondratyuka Str. 3 129515, Moscow Tel: (7-499) 553-6036

Address: Avenida Manoteras 38, Bloque B, Bajo B008. Madrid 28050, Spain Tel: (34-91)392-3754

Address: General Ali Rıza Gürcan Cad. Metropol Center No:31 A Blok K.13 D.52 Merter/Istanbul, Turkey Tel: (90-212) 482-0877

Address: 3 Percy Road, St. John's Park, Huntingdon, Cambs, Pe29

6SZ, United Kingdom Tel: (44-1480) 416-840

DISTRIBUTOR:

USAAddress: 800 MacArthur Blvd., Mahwah, NJ 07430-0619, USA
Tel: (1-800) 288 2121, (1-201) 995-8000





Mindray is listed on the NYSE under the symbol "MR" Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China Tel: +86 755 26582888 Fax: +86 755 26582680 E-mail: intl-market@mindray.com Website: www.mindray.com

